

Characterization of a New High Performance Inert Nebulizer for ICP Spectrometry

Jerry Dulude, Ryan Brennan, and Jol Desmarchelier

Glass Expansion

Characteristics of Concentric Glass Nebulizers

- High efficiency
- High precision
- Small droplet size
- Well understood
- Resistant to organic solvent
- Resistant to most acids
- Not resistant to HF

Nukiyama and Tanasawa equation*

$$d_{3,2} = \frac{585}{V} \left[\frac{\sigma}{\rho} \right]^{0.5} + 597 \left[\frac{\eta}{(\sigma\rho)^{0.5}} \right]^{0.45} \left[\frac{10^3 Q_l}{Q_g} \right]^{1.5}$$

$d_{3,2}$ = Sauter mean diameter - (μm)

V = Velocity difference of gas-liquid - (m/s)

σ = Surface tension - (dyn/cm)

ρ = Liquid density - (g/cm^3)

η = Liquid viscosity - (Poise or $\text{dyn}\cdot\text{s}/\text{cm}^2$)

Q_l = Volume flowrate, liquid - (cm^3/s)

Q_g = Volume flowrate, gas - (cm^3/s)

* S Nukiyama and Y Tanasawa, Trans Soc Mech Eng., (1938-1940)

Currently Available Concentric Inert Nebulizers

OpalMist

- PFA body and insert
- 3% physical reproducibility
- Up to 15% TDS
- Up to 75um particles
- Highest purity
- 0.05, 0.1, 0.2, 0.4, 0.6 & 2mL/min models



PolyCon

- Polyimide body and insert
- 2% physical reproducibility
- Up to 5% TDS
- Up to 75um particles (for high uptake model)
- 0.05, 0.1, 0.2, 0.4, 0.6, 0.8, 2 & 5mL/min



New DuraMist Nebulizer

- PEEK body and insert
- 2% physical reproducibility
- Up to 30% TDS
- Up to 75um particles (for larger uptake)
- 0.4 and 1.0mL/min uptake models
- Either 1.0 or 0.7 L/min argon flow



DuraMist Development Goals




- HF resistant concentric nebulizer
- Lower price than PolyCon or OpalMist
- Similar performance to PolyCon
- Tolerant of various sample types
- Approach performance specs of SeaSpray

Thermo ICAP 6300 Operating Conditions

RF Power	1350 W
Plasma gas	15 L/min
Aux gas	0.2 L/min
Neb gas	0.65 L/min
Replicates	3
Sample flush time	65 sec
Plasma view	auto
Max Integration	15 sec
Analysis pump rate	37 rpm
Pump tubing	Orange/white


IsoMist 21° C twister spray chamber, D-Torch.

Comparison of Precision




NCPN = non-concentric polymer nebulizer

Comparison of Intensity to SeaSpray



DuraMist Long-term Stability




ICP Figures of Merit

Abbreviation	Test	Diagnostic
Mg-Mg	Mg280(II)/Mg285(I)	Robustness
Cr-Cr	Cr267(II)/Cr357(I)	Atomization/ionization
Zn-Ba	Zn206(II)/Ba455(II)	Excitation
Mg	Mg285 (I) RSD	Neb. Efficiency
Ar	Ar404 (I) RSD	Stability
Zn	Zn206 (II) RSD	Stability
Ba	Ba455 (II) RSD	Stability


**Simple experiments for the control, the evaluation and the diagnosis of
Inductively coupled plasma sequential systems**

E. POUSSEL and J. M. MERMET
Spectrochimica Acta, 1993

Comparison of figures of merit (Higher is better)



Comparison of figures of merit (lower is better)




Comparison with other inert nebulizers

PE Optima 2100 Operating Conditions

RF Power	1400 W
Plasma gas	15 L/min
Aux gas	0.2 L/min
Neb gas	0.65 L/min
Read delay time	15 sec
Replicates	3
Min Integration	1 sec
Max Integration	10 sec
Source Equilibration delay	15 sec
Plasma View	Axial (Na_589 = radial)

IsoMist 21° C twister spray chamber, D-Torch


Comparison of Sensitivity




Tolerance to Dissolved Salts (20%)

- Natural aspiration of 20 % NaCl.
- Nebulizer argon flow monitored at constant pressure
- Total time of 20 minutes.


Tolerance to Dissolved Salts (20%)




DuraMist Reproducibility



Comparison of figures of merit (Higher is better)



Comparison of figures of merit (lower is better)



DuraMist Summary

- Intensities approach the best concentric glass nebulizer – the SeaSpray
- Excellent precision
- HF resistant
- Durable
- Cost effective