

Glass DC Nebulizers

SeaSpray

| | |
|------------------|--------------|
| 2mL/min uptake | A21-07-USS2 |
| 400µL/min uptake | A21-07-USS04 |

MicroMist

| | |
|------------------|--------------|
| 600µL/min uptake | A21-07-UM06 |
| 200µL/min uptake | A21-07-UM02 |
| 100µL/min uptake | A21-07-UM01 |
| 50µL/min uptake | A21-07-UM005 |



Conikal

| | |
|----------------|------------|
| 2mL/min uptake | A21-07-UC2 |
| 1mL/min uptake | A21-07-UC1 |

Slurry

| | |
|----------------|------------|
| 4mL/min uptake | A21-07-US6 |
|----------------|------------|

DC Nebulizers for HF

DuraMist

| | |
|------------------|-------------|
| 1mL/min uptake | A21-07-DM1 |
| 400µL/min uptake | A21-07-DM04 |



OpalMist

| | |
|------------------|---------------|
| 2mL/min uptake | A21-07-PFA2 |
| 600µL/min uptake | A21-07-PFA06 |
| 200µL/min uptake | A21-07-PFA02 |
| 100µL/min uptake | A21-07-PFA01 |
| 50µL/min uptake | A21-07-PFA005 |



VeeSpray

| | |
|----------------------|------------|
| 0.6 - 3mL/min uptake | A21-07-CV6 |
|----------------------|------------|



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

| | |
|----------------------|-------------|
| Glass spray chamber | KT-1012-XR |
| Quartz Spray Chamber | KT-1021Q-XR |
| PFA spray chamber | KT-1012P-XR |



Spray Chambers

Tracey

| | |
|---------------|---------------|
| 50mL Cyclonic | 20-809-0434HE |
|---------------|---------------|



Twister

| | |
|--|---------------|
| 50mL Cyclonic (Aqueous) With 7mm baffle | 20-809-0208HE |
|--|---------------|



Twister

| | |
|---|---------------|
| 50mL Cyclonic (Organics) With 4mm baffle | 20-809-0446HE |
|---|---------------|



Twister

| | |
|---------------------------|---------------|
| 50mL Cyclonic (Aux. Port) | 20-809-0004HE |
|---------------------------|---------------|



Cinnabar

| | |
|---------------|---------------|
| 20mL Cyclonic | 20-809-9121HE |
|---------------|---------------|



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

| | |
|---------------|---------|
| HydraMist Kit | KT-1162 |
|---------------|---------|

Spray Chambers for HF

Tracey TFE

| | |
|------------------------------------|-------------|
| 50mL Cyclonic | 20-809-2510 |
| Plugs directly into cassette mount | |



Torches

D-Torch 30-808-3300

Comprises:

Base and Inner Tube 31-808-3292
 Quartz Outer Tube 31-808-3294
 Retaining Ring 31-808-2818



D-Torch Options

Ceramic Outer Tube 31-808-2990
 Quartz Outer Tube (No Slot) 31-808-3407
 Copper Foils (PKT 12) for Ceramic Outer Tube 31-808-3356



Quartz Torch

Without slots 31-808-3254
 With single slot 31-808-3255
 With 3 slots 31-808-3256

Quartz Torch Options:

Copper Foil for Quartz Torch or D-Torch Quartz Outer Tube (pack of 12) 31-808-3289

Tapered Quartz Injectors

Quartz Injector 0.8mm 31-808-3257
 Quartz Injector 2.0mm 31-808-3258
 Quartz Injector 2.5mm 31-808-3730
 Quartz Injector 3.0mm 31-808-3259
 Alumina Injector 1.2mm 31-808-3287
 Alumina Injector 2.0mm 31-808-3286



Sapphire Injectors

Sapphire Injector 2mm 31-808-3528

Purge Windows

Axial 17mm disk 31-808-0551
 Axial 32mm disk 31-808-2649
 Radial (Post Nov. 2004) 31-808-2788

Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1269



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Fitting only connects to 1.3mm and 1.6mm Autosampler Tubing



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF solutions 60-808-1150



Autosampler Probes

(See page 14 for complete list of probes)

Suitable for Cetac ASX-200/500/800 Series

Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793

Suitable for PE S10 or AS93+

PTFE Encapsulated Carbon Fibre Probe 0.25mm ID 70-803-1071

PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-1440

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0816

