

Glass DC Nebulizers

Conikal
1mL/min uptake A41-1-UC1S

SeaSpray
2mL/min uptake A41-1-USS2S

Slurry
4mL/min uptake A41-1-US6

MicroMist
600µL/min uptake A41-07-UM06S
200µL/min uptake A41-07-UM02
100µL/min uptake A41-07-UM01
50µL/min uptake A41-07-UM005



DC Nebulizers for HF

OpalMist
1mL/min uptake A41-07-PFA1S
600µL/min uptake A41-07-PFA06S
200µL/min uptake A41-07-PFA02
100µL/min uptake A41-07-PFA01
50µL/min uptake A41-07-PFA005



DuraMist
1mL/min uptake A41-07-DM1S
400µL/min uptake A41-07-DM04



VeeSpray
0.6 - 3mL/min uptake A41-1-CV6



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

Glass spray chamber KT-1067-XR
Quartz spray chamber KT-1067Q-XR
PFA spray chamber KT-1067P-XR



Spray Chambers

Tracey
50mL Cyclonic 20-808-8881HE
With B14 gravity drain



Twister
50mL Cyclonic 20-809-9119HE
With B14 gravity drain



Twister
50mL Cyclonic 20-809-0222HE



Cinnabar
20mL Cyclonic 20-809-0164HE



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Hydride Kit KT-1157



Drain Device 21-808-8214



Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

Glass Transfer Tube 86mm 21-809-2846



Spray Chambers for HF

Tracey TFE

50mL Cyclonic 20-809-2508
 With B14 gravity drain
 Suitable for use with HF



HF Drain Device

21-809-2540



Mounting bracket for HF Spray Chamber

21-809-2839



Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

PFA Transfer Tube 82.5mm 31-808-3071



Torches

Semi Demountable Torch (HF)

30-808-8039

Suitable for ICPS-7500/8100 and ICPE-9000/9800

Comprises:

Torch Body 31-808-1258
 Torch Adaptor 31-808-1260
 Tapered Alumina Injector 1.5mm 31-808-0013
 Joint Clip 14/23 KCP-14
 GazFit Connectors (PKT 4) for 4mm OD side arms GAZ-04



Injectors for Semi Demountable Torch

Alumina Injector 1.0mm 31-808-0835
 Alumina Injector 1.5mm 31-808-0013
 Alumina Injector 1.8mm 31-808-0015
 Alumina Injector 2.0mm 31-808-0027
 Alumina Injector 2.4mm 31-808-0385
 Sapphire Injector 2.0mm 31-808-3001



Torches

D-Torch

30-808-3910

Suitable for ICPE-9000/9800

Comprises:

Torch Body 31-808-3998
 Torch Adaptor 31-808-3072
 Injector Ferrule 6.0mm 31-808-2962
 Quartz Outer Tube 31-808-3912
 Outer Tube Retainer 31-808-2818
 GazFit Connectors (PKT 4) for 4mm OD side arms GAZ-04

D-Torch Options:

Quartz Outer Tube 31-808-3912


Tapered Injectors for D-Torch

Alumina Injector 1.5mm 31-808-3964
 Quartz Injector 1.0mm 31-808-3965
 Quartz Injector 1.5mm 31-808-3966
 Quartz Injector 1.8mm 31-808-3967



Quartz Torches

Quartz Torch Standard 1.5mm Injector (Suitable for 7500/8100) 30-808-1179
 Quartz Torch for Organics, 1.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0519
 Quartz Torch for High TDS, 1.8mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0518
 Quartz Torch for High TDS, 2.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0564
 Quartz Mini Torch 1.5mm Injector (Suitable for ICPE-9000/9800) 30-807-0534



RF Coils

Copper	70-900-6100C
Silver	70-900-6100S
Gold	70-900-6100G



RF Coil Installation Kit 70-900-6100T

Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1312



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Autosampler Probes

(See page 14 for complete list of probes)

For Shimadzu® ASC-6100F/9800

PTFE Sheathed Carbon Fibre Probe
0.3mm ID with UniFit 70-803-1200

For Shimadzu® AS10

PTFE Coated Vitcon O-rings
(PKT 10) With UniFit 70-803-1477

