

## Glass DC Nebulizers

### SeaSpray

400µL/min uptake      A13-1-USS04

### MicroMist

400µL/min uptake      A13-1-UM04

200µL/min uptake      A13-1-UM02

100µL/min uptake      A13-1-UM01

50µL/min uptake      A13-1-UM005

### Conikal

1mL/min uptake      A13-1-UC1



## Quartz DC Nebulizers

### SeaSpray

1mL/min uptake      A13-1-QSS1



## DC Nebulizers for HF

### OpalMist

400µL/min uptake      A13-1-PFA04

200µL/min uptake      A13-1-PFA02

100µL/min uptake      A13-1-PFA01

50µL/min uptake      A13-1-PFA005



## DC Nebulizers for HE-SIS

### MicroMist HE

200µL/min uptake      C13-1-UFT02

200µL/min uptake      E13-1-UFT02



## Peltier Controlled Cyclonic Spray Chambers

Glass Spray Chamber      KT-1116

Quartz Spray Chamber      KT-1116Q

PFA Spray Chamber      KT-1116P



## Spray Chambers

### High Efficiency Sample Introduction System

HE-SIS for Agilent ICP-MS      KT-1155



### \*Twister

50mL Cyclonic with Aux.      20-809-2574HE

\*Requires adaptor KT-1104.



### Quartz Scott Style

For 7900 UHMI      20-809-4144



## Spray Chambers for HF

### \*Tracey PFA44

44ml Cyclonic with Aux.      20-809-2578

Spray Chamber Adaptor      KT-1104

\*Requires adaptor KT-1104.



## Torches

### Quartz Torches

with 1.0mm Injector      30-807-0558

with 1.5mm Injector      30-807-0557

with 2.0mm Injector      30-807-3608

with 2.5mm Injector      30-807-0556



### Torch Options

Quartz Torch Bonnet      31-808-2577

Quartz Connector with Dilution Port      31-808-3217

Quartz Connector with 45deg Dilution Port      31-808-3280

Straight Quartz Connector      31-808-3288

Clamp (PKT 2)      70-803-0846



**Torches**

(continued)

**Semi-Demountable Torch**

with 2.4mm Quartz Injector 30-808-3853



**Torch Options**

Semi-Demountable Torch body 31-808-3857



1.0mm Tapered Quartz Injector 31-808-4172

1.5mm Tapered Quartz Injector 31-808-4173



2.4mm Tapered Quartz Injector 31-808-4168

1.0 mm Tapered Alumina Injector 31-808-4169



1.5 mm Tapered Alumina Injector 31-808-4170

2.4mm Tapered Alumina Injector 31-808-3855

2.4mm Sapphire Capillary Injector 31-808-4159



Injector Adaptor 31-808-4129

Viton O-rings for FDT/SDT ball joints (PKT 10) 70-V-011



**RF Coils**

Copper 70-900-7800C



Silver 70-900-7800S



Gold 70-900-7800G




**RF Coil Installation Kit** 70-900-7500T



**Cones**

**Sampler Cones**

Nickel Sampler AT7701-Ni



Nickel Sampler with Nickel plated base AT7701-Ni/Ni



Aluminium Sampler AT7701-Al



Platinum Sampler AT7706-Pt

Platinum Sampler (18 mm Insert) AT7706A-Pt




Platinum Sampler with Nickel plated base AT7706-Pt/NiP



Platinum Sampler with solid Nickel base AT7706-Pt/Ni

Platinum Sampler (15mm insert) with Nickel base AT7706B-Pt/Ni



Platinum Sampler with Nickel plated base (Extended Life) AT7706L-Pt/NiP

**Skimmer Cones**

Nickel Skimmer for x-lens AT7902X-Ni



Aluminium Skimmer for x-lens AT7902X-Al



Platinum Skimmer for x-lens AT7908X-Pt



## Cones

(continued)

### Skimmer Cones (continued)

Platinum Skimmer with Nickel base for x-lens AT7908X-Pt/Ni



Nickel Skimmer for s-lens AT7702S-Ni

Copper Skimmer for s-lens AT7702S-Cu



Platinum Skimmer for s-lens AT7708S-Pt

Platinum Skimmer with Nickel base for s-lens AT7708S-Pt/Ni



Platinum Skimmer with Nickel base for m-lens AT8908M-Pt/Ni

### Cone Options

Skimmer Cone Guard 70-803-1004

Long Life Shield Plate AT5004

Graphite Gasket (PKT 3) AT7703

Retaining Ring AT7704




Brass Skimmer Adaptor base for m lens AT7905M-BR

Brass Skimmer Adaptor base for x lens AT7905X-BR

Stainless Steel Skimmer Adaptor base for x lens AT7905X-SS

Captive Screw M4 for Skimmer Base (PKT 3) 70-803-2061



## Common Accessories

### Elegra Argon Humidifier

Elegra Dual Argon Humidifier 70-803-1273



Elegra Stand 70-803-1581



## Common Accessories

(continued)

### Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



### Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



## Autosampler Probes (See page 20-21 for complete list of probes)

### For Agilent I-AS

PTFE Encapsulated Carbon Fibre Probe 0.25mm ID 60-703-1010

PTFE Encapsulated Carbon Fibre Probe 0.3mm ID 60-703-1009



### For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793



### For SPS 3/SPS 4

Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-1957

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-0908

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0853

