

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

Product List for Spectro (Ametek) ICP-OES: Spectroblue EOP and TI



GLASS EXPANSION
Quality By Design

INTERNATIONAL

15 Batman Street
West Melbourne
Vic 3003, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

Spectro (Ametek) ICP-OES: Spectroblue EOP and TI

Glass Nebulizers

Conikal
2mL/min uptake [ARG-07-UC2](#)

SeaSpray
2mL/min uptake [ARG-07-USS2](#)

Slurry
4mL/min uptake [ARG-07-US6](#)

MicroMist
600µL/min uptake [ARG-07-UM06](#)
200µL/min uptake [ARG-07-UM02](#)
100µL/min uptake [ARG-07-UM01](#)
50µL/min uptake [ARG-07-UM005](#)



Nebulizers suitable for use with HF

OpalMist
2mL/min uptake [ARG-07-PFA2](#)
600µL/min uptake [ARG-07-PFA06](#)
200µL/min uptake [ARG-07-PFA02](#)
100µL/min uptake [ARG-07-PFA01](#)
50µL/min uptake [ARG-07-PFA005](#)



DuraMist
1mL/min uptake [ARG-07-DM1](#)
400µL/min uptake [ARG-07-DM04](#)



VeeSpray - Ceramic
Standard [ARG-07-CV6](#)



IsoMist


IsoMist Programmable Temperature Spray Chamber

With glass spray chamber [KT-1027](#)
With quartz spray chamber [KT-1027Q](#)
With PFA spray chamber [KT-1027P](#)



Spray Chambers

Twister
50mL Cyclonic [20-809-2951HE](#)



Quartz Adaptor
[21-809-3772](#)



Quartz Adaptor with Aux Port
[21-809-3771](#)



Tracey TFE
50mL Cyclonic [20-809-3717](#)
Suitable for use with HF.
Requires Transfer Tube [21-809-3726](#).



PuraFlex Transfer Tube
[21-809-3726](#)



Torches

Quartz Torch
[31-808-3394](#)



Injectors
Quartz Injector 1.2mm [31-808-3410](#)
Quartz Injector 1.8mm [31-808-3391](#)
Quartz Injector 2.0mm [31-808-3411](#)
Quartz Injector 3.0mm [31-808-3412](#)



Alumina Injector 1.8mm [31-808-3409](#)
Alumina Injector 2.5mm [31-808-3390](#)

