

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

Product List for Thermo ICP-OES: iCAP 6000/7000 Duo



GLASS EXPANSION
Quality By Design

INTERNATIONAL

15 Batman Street
West Melbourne
Vic 3003, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

Glass Nebulizers

Conikal	
2mL/min uptake	ARG-07-UC2
SeaSpray	
2mL/min uptake	ARG-07-USS2
Slurry	
4mL/min uptake	ARG-07-US6
MicroMist	
600µL/min uptake	ARG-07-UM06
200µL/min uptake	ARG-07-UM02
100µL/min uptake	ARG-07-UM01
50µL/min uptake	ARG-07-UM005



Nebulizers suitable for use with HF

OpalMist	
2mL/min uptake	ARG-07-PFA2
600µL/min uptake	ARG-07-PFA06
200µL/min uptake	ARG-07-PFA02
100µL/min uptake	ARG-07-PFA01
50µL/min uptake	ARG-07-PFA005
DuraMist	
1mL/min uptake	ARG-07-DM1
400µL/min uptake	ARG-07-DM04
VeeSpray - Ceramic	
Standard	ARG-07-CV6



IsoMist



IsoMist Programmable Temperature Spray Chamber for Thermo iCAP 6000/7000 Duo

With glass spray chamber	KT-1015
With quartz spray chamber	KT-1015Q
With PFA spray chamber	KT-1015P



Spray Chambers

Tracey	
50mL Cyclonic	20-809-2665HE
Twister	
50mL Cyclonic	20-809-2951HE
Cinnabar	
20mL Cyclonic	20-809-3220HE



Tracey TFE	
50mL Cyclonic (Suitable for HF)	20-809-3469



Adaptor	
Right Angle With Ball Joint	21-809-2761



Torches

D-Torch Demountable Torch

D-Torch with Quartz Outer Tube	30-808-2844
--------------------------------	-----------------------------

Comprises:	
Base and Inner Tube	31-808-2932
Quartz Outer Tube	31-808-2840
Retaining Ring	31-808-2707
Injector Adaptor	31-808-2895
Injector Ferrule	31-808-2897



Torches (continued)

Injectors for D-Torch & EMT Torch

Tapered Quartz Injector 1.0mm [31-808-2833](#)

Tapered Quartz Injector 1.5mm [31-808-2834](#)

Capillary Quartz Injector 2.0mm [31-808-2836](#)

Capillary Quartz Injector 2.5mm [31-808-3246](#)

Capillary Quartz Injector 3.0mm [31-808-3247](#)

Alumina Injector 1.0mm [31-808-3279](#)

Alumina Injector 2.0mm [31-808-2835](#)

Sapphire Injector 2.0mm [31-808-2977](#)



Ceramic Outer Tube and Optic Fibre for D-Torch

[31-808-3079](#)

Comprises:

Ceramic Outer Tube [31-808-2863](#)

Optic Fibre 170mm [31-808-2933](#)

Optic Fibre 450mm [31-808-3117](#)



Quartz Torch

EMT Design [31-808-2811](#)



Quartz Torch

Old Design [31-808-2676](#)



Injectors for Old Design Torch

Tapered Quartz Injector 1.0mm [31-808-2719](#)

Tapered Quartz Injector 1.5mm [31-808-2702](#)

Capillary Quartz Injector 2.0mm [31-808-2675](#)

Tapered Alumina Injector 1.5mm [31-808-2787](#)

Alumina Injector 2.0mm [31-808-2786](#)

Sapphire Injector 2.0mm [31-808-2804](#)



Other Products

Duo Periscope Window [31-808-2677](#)



Quartz Window [70-803-1116](#)



RF Coils

RF Coil

Gold, Teflon Coated [70-900-4002G](#)



RF Coil

Silver, Teflon Coated [70-900-4002S](#)



Installation Kit

[70-900-4002T](#)

