

# ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

## Product List for Thermo ICP-OES: iCAP 6000/7000 Radial



**GLASS EXPANSION**  
Quality By Design

### INTERNATIONAL

15 Batman Street  
West Melbourne  
Vic 3003, Australia

Telephone: +61 3 9320 1111  
Facsimile: +61 3 9320 1112  
Email: [enquiries@geicp.com](mailto:enquiries@geicp.com)

### AMERICAS

4 Barlows Landing Road, Unit 2A  
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097  
Telephone: 508 563 1800  
Facsimile: 508 563 1802  
Email: [geusa@geicp.com](mailto:geusa@geicp.com)

## Glass Nebulizers

**Conikal**  
2mL/min uptake [ARG-07-UC2](#)

**SeaSpray**  
2mL/min uptake [ARG-07-USS2](#)

**Slurry**  
4mL/min uptake [ARG-07-US6](#)

**MicroMist**  
600µL/min uptake [ARG-07-UM06](#)

200µL/min uptake [ARG-07-UM02](#)

100µL/min uptake [ARG-07-UM01](#)

50µL/min uptake [ARG-07-UM005](#)



## Nebulizers suitable for use with HF

**OpalMist**  
2mL/min uptake [ARG-07-PFA2](#)

600µL/min uptake [ARG-07-PFA06](#)

200µL/min uptake [ARG-07-PFA02](#)

100µL/min uptake [ARG-07-PFA01](#)

50µL/min uptake [ARG-07-PFA005](#)



**DuraMist**  
1mL/min uptake [ARG-07-DM1](#)

400µL/min uptake [ARG-07-DM04](#)



**VeeSpray - Ceramic**  
Standard [ARG-07-CV6](#)



## IsoMist

### IsoMist Programmable Temperature Spray Chamber for Thermo iCAP 6000/7000 Radial

With glass spray chamber [KT-1021](#)

With quartz spray chamber [KT-1021Q](#)

With PFA spray chamber [KT-1021P](#)




## Spray Chambers

**Tracey**  
50mL Cyclonic [20-809-2665HE](#)



**Twister**  
50mL Cyclonic [20-809-2951HE](#)



**Tracey TFE**  
50mL Cyclonic [20-809-9198](#)  
Suitable for use with HF.



**Straight Ball Joint Adaptor**  
[21-809-2723](#)



## Torches

**D-Torch Demountable Torch**  
[30-808-2877](#)

**Comprises:**  
Base and Inner Tube [31-808-2932](#)  
Quartz Outer Tube [31-808-2878](#)  
Retaining Ring [31-808-2707](#)  
Injector Adaptor [31-808-2895](#)  
Injector Ferrule [31-808-2897](#)



**Options:**  
Ceramic Outer Tube [31-808-2881](#)



### Injectors for D-Torch & EMT Torch


Tapered Quartz Injector  
1.0mm (EMT) [31-808-2833](#)

Tapered Quartz Injector  
1.5mm (EMT) [31-808-2834](#)

Capillary Quartz Injector  
2.0mm (EMT) [31-808-2836](#)

Alumina Injector 2.0mm (EMT) [31-808-2835](#)

Sapphire Injector 2.0mm [31-808-2977](#)



# Thermo ICP-OES: iCAP 6000/7000 Radial

## Torches (continued)

### Quartz Torches

Quartz Torch (EMT) [31-808-2810](#)



Quartz Torch (Old Design) [31-808-2692](#)



### Injectors for Old Design Semi Demountable Torch

Tapered Quartz Injector 1.0mm [31-808-2719](#)



Tapered Quartz Injector 1.5mm [31-808-2702](#)

Capillary Quartz Injector 2.0mm [31-808-2675](#)

Tapered Alumina Injector 1.5mm [31-808-2787](#)

Alumina Injector 2.0mm [31-808-2786](#)

Sapphire Injector 2.0mm [31-808-2804](#)



## RF Coils

### RF Coil

Gold, Teflon Coated [70-900-4002G](#)



### RF Coil

Silver, Teflon Coated [70-900-4002S](#)



### Installation Kit

[70-900-4002T](#)

